

## Determining the gravitational acceleration with a reversible pendulum

### Objects of the experiments

- Measuring the oscillation periods  $T_1$  and  $T_2$  of a reversible pendulum for two suspension points.
- Tuning the reversible pendulum to the same oscillation period.
- Determining the gravitational acceleration from the oscillation period and the reduced length of pendulum.




Fig. 1 Oscillations of a reversible pendulum around the suspension points  $H_1$  and  $H_2$ .

### Principles

#### Compound pendulum:

If a compound pendulum oscillates around its rest position with small deflections  $\varphi$ , the equation of motion is:

$$J \cdot \ddot{\varphi} + m \cdot s \cdot g \cdot \varphi = 0 \quad (I)$$

$J$ : moment of inertia around the axis of oscillation,  
 $s$ : distance between the axis of oscillation and the centre of mass,  
 $g$ : gravitational acceleration,  $m$ : mass of the pendulum

The reduced length of the compound pendulum is defined as the quantity

$$s_r = \frac{J}{m \cdot s} \quad (II)$$

because its oscillation period

$$T = 2\pi \sqrt{\frac{s_r}{g}} \quad (III)$$

corresponds to that of a simple pendulum with the length  $s_r$ .

The moment of inertia  $J$  of the compound pendulum is, according to the parallel axis theorem,

$$J = J_S + m \cdot s^2 \quad (IV)$$

$J_S$ : moment of inertia around the centre of mass axis

Therefore the reduced length of pendulum is

$$s_r = \frac{J_S}{m \cdot s} + s \quad (V)$$

#### Reversible pendulum:

The reversible pendulum is a particular type of the compound pendulum. There are two edges  $H_1$  and  $H_2$  that allow to choose the suspension point. Two masses  $m_1 = 1000$  g and  $m_2 = 1400$  g on the straight line  $H_1H_2$  can be shifted so that the oscillation period is tunable. The goal of the tuning is to achieve equal oscillation periods around both edges. In this case, the reduced length of pendulum is equal to the distance  $d = 99.4$  cm between the edges. This latter statement can be understood from the following consideration:

**Apparatus**

1 reversible pendulum . . . . .	346 111
1 steel tape measure, 2m . . . . .	311 77
1 stopclock I, 30s/15 min . . . . .	313 07

Because of Eq. (II) the pendulum oscillates around both edges  $H_1$  and  $H_2$  at the same oscillation period  $T_1 = T_2$  if the reduced lengths of pendulum  $s_{r,1}$  and  $s_{r,2}$  agree. In this case we have

$$\frac{J_S}{m \cdot s_1} + s_1 = \frac{J_S}{m \cdot s_2} + s_2 \quad (VI).$$

The relation between the distances  $s_1$  and  $s_2$  between the two edges and the centre of mass is given by

$$s_1 + s_2 = d \quad (VII)$$

because the centre of mass  $S$  is located on the connecting line  $H_1H_2$  for reasons of symmetry.

From Eqs. (VI) and (VII) a quadratic equation is obtained for determining  $s_1$ . Its solution is

$$s_1 = \frac{d}{2} \pm \sqrt{\frac{d^2}{4} - \frac{J_S}{m}} \quad (VIII).$$

If  $s_1$  is inserted in Eq. (V) instead of  $s$ , the reduced length of pendulum of the tuned reversible pendulum is obtained:

$$s_r = d \quad (IX).$$

**Determining the gravitational acceleration:**

For the tuned reversible pendulum

$$T^2 = 4 \cdot \pi^2 \cdot \frac{d}{g} \quad (X)$$

follows from Eqs. (III) and (IX). The distance  $d$  between the edges is known to a high precision. After the oscillation period  $T$  of the tuned pendulum has been determined, the gravitational acceleration  $g$  is obtained from Eq. (X).

**Setup**

The experimental setup is illustrated in Fig. 2.

- If necessary, assemble the reversible pendulum according to the instruction sheet, and mount the wall holder on a stable, non-vibrating wall.
- Make marks on the pendulum rod at  $x_1 = 25$  cm and  $x_2 = 10$  cm, 15 cm, 20 cm etc. (starting from the edge  $H_1$ , take into account the radius  $r = 5$  cm of the masses, see Fig. 3)
- Fix the mass  $m_1$  at the position  $x_1 = 25$  cm.





Fig. 2 Experimental setup for determining the gravitational acceleration with a reversible pendulum.

Fig. 3 Distance marks  $x_1$  and  $x_2$  on the pendulum rod.



**Carrying out the experiment**

- Carefully hang the reversible pendulum on the edge bearing of the wall holder with the edge  $H_1$  and fix the mass  $m_2$  at the distance  $x_2 = 50$  cm.
- Cautiously deflect the lower end of the reversible pendulum parallel to the wall, and make it oscillate, if possible, without any vibrations.
- Measure the length of 50 oscillations and take it down as measured value  $50 \cdot T_1$ .

- Hang the reversible pendulum on the edge bearing with the edge  $H_2$ , and measure the period  $50 \cdot T_2$ .
- Slide the mass  $m_2$  to the position  $x_2 = 55$  cm, and measure  $50 \cdot T_2$  at first and then  $50 \cdot T_1$ .
- Slide the mass  $m_2$  towards the edge  $H_2$  in steps of 5 cm; each time measure the two oscillation periods. Plot  $T_1^2$  and  $T_2^2$  as functions of  $x_2$  and, if necessary, repeat the measurement of the oscillation periods.
- Next slide the mass  $m_2$  towards the edge  $H_1$  in steps of 5 cm starting from  $x_2 = 45$  cm, and measure the two oscillation periods each time.

Evaluation and results

The two measured curves  $T_1^2(x_2)$  and  $T_2^2(x_2)$  intersect near  $x_2 = 30$  cm and  $x_2 = 65$  cm (see Fig. 4). The enlarged sections in Figs. 5 and 6 shown that the curves intersect at  $T^2 = 4.039$  s<sup>2</sup> and  $T^2 = 4.014$  s<sup>2</sup>, respectively. With the mean value  $T^2 = 4.027$  s<sup>2</sup> Eq. (X) gives

$$g = \frac{4 \cdot \pi^2 \cdot d}{T^2} = 9.74 \frac{\text{m}}{\text{s}^2}$$

Measuring example

Table 1: Oscillation periods  $T_1$  and  $T_2$  around the edges  $H_1$  and  $H_2$ , respectively, as functions of the distance  $x_2$  between the mass  $m_2$  and the edge  $H_1$ .

$\frac{x_2}{\text{cm}}$	$\frac{50 \cdot T_1}{\text{s}}$	$\frac{T_1^2}{\text{s}^2}$	$\frac{50 \cdot T_2}{\text{s}}$	$\frac{T_2^2}{\text{s}^2}$
20	106.0	4.494	101.9	4.153
25	103.1	4.252	101.2	4.097
30	100.8	4.064	100.6	4.048
35	99.4	3.952	100.1	4.008
40	98.8	3.905	99.8	3.984
45	98.2	3.857	99.6	3.968
50	97.9	3.834	99.5	3.960
55	98.3	3.865	99.6	3.968
60	99.1	3.928	99.8	3.984
65	99.9	3.992	100.0	4.000
70	101.1	4.088	100.7	4.056
75	102.2	4.178	101.7	4.137
80	103.2	4.260	102.2	4.178
85	104.8	4.393	103.6	4.293




Fig. 5 Enlarged section of Fig. 4 around  $x_2 = 30$  cm with a non-linear interpolation of the measured values. Point of intersection: ( $x_2 = 31$  cm,  $T^2 = 4.039$  s<sup>2</sup>).

Fig. 4 Squared oscillation periods around the edges  $H_1$  (●) and  $H_2$  (■) as functions of the distance  $x_2$  between the mass  $m_2$  and the edge  $H_1$ .




Fig. 6 Enlarged section of Fig. 4 around  $x_2 = 65$  cm with a non-linear interpolation of the measured values. Point of intersection: ( $x_2 = 66$  cm,  $T^2 = 4.014$  s<sup>2</sup>).

