

## Determination of the standard electrode potentials of non-metals (Using the electrochemistry demonstration unit)

### Aims of the experiment

- To produce a standard hydrogen electrode.
- To recognise that even non-metals form an electrode potential.
- To measure the standard electrode potentials of halogens as an example on non-metals.
- To carry out the electrolysis of various solutions.
- To apply redox reactions.

### Principles

Electrochemistry deals with chemical reactions in which electrical current is generated or must be supplied. These are mostly redox reactions in which one reaction partner releases electrons and is oxidised. The other reaction partner accepts electrons and is reduced. The name redox results from the fact that both individual reactions always occur simultaneously.

Not all substances have the same tendency to release or to accept electrons. They can therefore be arranged as a series in order of voltage. The redox pairs are arranged in order of their individual standard electrode potential.

To determine the electrode potentials of corresponding redox pairs, half-cells are considered. Non-metals themselves are not electrically conductive and are often not able


to form an electrode as a solid substance. This is particularly obvious with gaseous non-metals, such as the halogens. For this reason, in this experiment a graphite electrode will be immersed in the individual halogen salt solutions. The graphite electrode is electrically conductive but chemically inert. Alternatively an inert metal electrode, such as a platinum electrode, could be used. This is immersed in the relevant salt solution and is surrounded by the gas of the non-metal to be determined.

Electrode potentials can be compared with one another if they are measured with reference to the standard hydrogen electrode. This consists of a platinum sheet which is surrounded by hydrogen gas. The solution that the platinum sheet is immersed in is a 1 M hydrochloric acid solution. If standard values of the electrode potentials are to be measured, this must be done under standard conditions, i.e. at 25



Fig. 1: Set-up of the experiment.

°C and 1.013 bar. The standard electrode potential of the hydrogen electrode was set arbitrarily to 0 V in 1912. The electrode potential of the standard hydrogen electrode is based on the following reaction:



The sheet platinum is used because hydrogen also cannot form an electrode and does not conduct electrically. An acid serves as the electrolyte solution, so that an equilibrium can be formed between the  $\text{H}^+$  ions of the acid and the hydrogen adsorbed on the platinum.

By definition, the standard electrode potentials of substances which release electrons to the standard hydrogen electrode are given a negative sign. In the case of substances which take on electrons from the standard hydrogen electrode, the standard electrode potentials are given a positive sign.

If one now compares the potential differences of various non-metals with the standard hydrogen electrode, an electrochemical series can be set up based on the standard potentials determined. The electrochemical series for non-metals with reference to the standard hydrogen electrode can be integrated into the electrochemical series for metals (see Fig. 2).




Fig. 2: Standard electrode potentials of metals and non-metals compared with the standard hydrogen electrode with the potential  $E^0 = 0.00 \text{ V}$ .

With the help of this electrochemical series, redox reactions can be predicted and quantitative conclusions can be drawn about them. For example, it can be predicted whether a reaction will take place voluntarily or what voltage will be needed to force a reaction to take place.

In this experiment, the standard electrode potentials of the non-metals chlorine, iodine and bromine will be determined. To achieve this, they will be connected in sequence with a standard hydrogen electrode to form an electrochemical cell. All gases will be produced by electrolysis at the electrodes at the start of the experiment. This is sufficient to enable the electrode potentials to be read on the electrochemistry demonstration unit.

### Risk assessment

Avoid contact of hexachloroplatinic acid with the skin! Hexachloroplatinic acid must on no account be allowed to enter the drains as it is highly toxic to the environment.

During the electrolysis of chloride, bromide and iodide solutions, the elemental halogens chlorine, bromine and iodine are produced in small amounts. Therefore only carry out the experiment in a well-ventilated place.

### Hexachloroplatinic acid



Signal word:  
Hazard

#### Hazard statements

H301 Toxic if swallowed.  
H314 Causes severe skin burns and eye damage.  
H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
H317 May cause allergic skin reactions.

#### Precautionary statements

P260 Do not inhale dust/fume/gas/mist/vapours/spray.

P301+P310 IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.

P303+P361+P353 IF ON SKIN (or hair): Remove/take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338 If in eyes: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

P405 Store locked up.

P501 Dispose of contents/container according to local/regional/national/international regulations.

### Equipment and chemicals

1	Electrochemistry demonstration unit, CPS.664 4071
1	Electrochemistry accessories set.....664 401
1	Panel frame C50, two-level, for CPS .....666 425
1	Table for electrochemistry, CPS .....666 472
4	Beaker, DURAN, 150 mL, tall .....602 032
1	Compact balance 200 g: 0.01 g .....667 7977
1	Measuring cylinder 100 mL, w. plastic base.665 754
1	Hydrochloric acid 1 mol/L, 1 L.....674 6910
1	Sodium chloride, 250 g .....673 5700
1	Potassium bromide, 250 g .....672 4920
1	Potassium iodide, 50 g.....672 6620
1	Hexachloroplatinic acid, 5 g .....672 1901
1	Water, pure, 1L .....675 3400

### Set-up and preparation of the experiment

#### Set-up of the apparatus

1. Insert the electrochemistry demonstration unit into the upper panel frame and supply with electrical power.
2. Place the table for electrochemistry, CPS into the frame below (see Fig. 1).






electrolysis. This manifests itself in the formation of small bubbles at the platinum electrode. After removing the cables from the power supply (5), the value displayed on the universal measuring instrument remains constant after a brief waiting period.

Reactions are also visible on the side of the non-metal half-cells. The smell of chlorine is noticeable during the electrolysis of the chloride solution. In the half-cell with the bromide solution, a brownish colour can be observed. In the half-cell with the iodide solution, a yellowish colour can be seen. Here we are dealing with the elemental forms of chlorine, bromine and iodine.

## Results

The following reactions take place:



Here, the equilibrium lies on the side of the molecular compounds for all non-metals.

The standard electrode potentials obtained in this experiment with the help of a platinised platinum electrode as a hydrogen electrode are as follows:

$$E^\circ(\text{Cl}_2/2 \text{Cl}^-) = + 1.32 \text{ V (literature value: + 1.36 V)}$$

$$E^\circ(\text{Br}_2/2 \text{Br}^-) = + 1.03 \text{ V (literature value: + 1.07 V)}$$

$$E^\circ(\text{I}_2/2 \text{I}^-) = + 0.52 \text{ V (literature value: + 0.58 V)}$$

The values measured scarcely differ from the theoretical values found in the literature.

The electrode potentials of the halogens fall with increasing atomic number. The higher the atomic number of the halogen, the lower is its electrode potential, and therefore also its power of oxidation.

## Cleaning and disposal

The solutions of the non-metals cannot be reused, as the solutions in the two cell troughs begin to mix with one another during the experiment. Place the remaining solutions of potassium bromide and potassium iodide into a collection container designated for inorganic salts waste for disposal. Sodium chloride solution can be disposed of in the laboratory drain.

The platinising solution can be stored in a labelled container and used again.