

Common properties of bodies

Density

Determining the density of air

Objects of the experiment

1. Determining the mass and volume of air
2. Calculating the density of air

Setup



Measuring example

| | |
|--|--------|
| Mass m_1 in g | 290.3 |
| Mass m_2 in g | 289.1 |
| Mass m_L in g | 1.2 |
| Volume V_L in cm^3 | 1000 |
| Density ρ_L in g/cm^3 | 0.0012 |

Evaluation

Like the density of solid bodies and liquids, the density of air can be determined from the mass and the volume:

$$\rho = \frac{m}{V}.$$

The density of air calculated from the measured values is: $\rho_L = 0.0012 \text{ g}/\text{cm}^3$.

The table value is: $\rho_L = 0.00129 \text{ g}/\text{cm}^3$.

Apparatus

| | |
|--|----------|
| 1 sphere with 2 stopcocks, glass, 1 l | 379 07 |
| 1 Support ring for round flask, 250 ml, cork | 667 072 |
| 1 Hand vacuum pump..... | 375 58 |
| 1 Single-pan suspension balance 610 Tara..... | 315 23 |
| 1 Measuring beaker, 1000 ml | 604 211 |
| 1 Measuring cylinder, 1000 ml, with plastic base. | 665 757 |
| 1 Silicone tubing, 6 mm/ diam., 1 m..... | 604 432 |
| 1 Water, pure, 5 l | 675 3410 |

Carrying out the experiment

- Put the sphere with the support ring on the pan.
- Determine the mass m_1 of the sphere with support ring.
- Connect the vacuum pump to a hose nipple of the sphere, and open the stopcock on this side.
- Close the stopcock on the other side of the sphere.
- Evacuate the sphere, close the stopcock, and disconnect the vacuum pump from the sphere.
- Put the evacuated sphere with support ring on the pan, and determine the mass m_2 .
- Calculate the mass of the air by forming the difference of m_1 and m_2 .
- Determine the volume of the air the same way as in D 1.1.1.6.
- Calculate the density of air.