Hooke's law - Stand setup

### Object of the experiment

1. Investigating the relation between the force acting on a helical spring and the spring elongation

# Setup



## **Apparatus**

1 helical spring, 32 N/m	352 12
1 Weight, 0.1 kg	683 10
1 Weight, 0.2 kg	683 11
1 Weight, 0.5 kg	315 38
1 Weight, 1 kg	315 39
1 Metal rule 1 m	311 02
1 Pointer, pair	301 29
1 Stand base, V-shape, large	300 01
1 Stand rod, 100 cm, 12 mm diam	300 44
1 Clamp with hook	301 08

### Carrying out the experiment

- Suspend the helical spring from the clamp with hook.
- Mark the lower edge of the unloaded helical spring with a pointer and define this point to be the zero.
- Suspend the weights 0.1 kg and 0.2 kg from the helical spring as acting force *F*.
- Mark the lower edge of the helical spring, which now is elongated, with the second pointer.


- Measure the spring elongation s between the two pointers.
- Enhance the acting force *F* by suspending further weights, and measure the corresponding spring elongations.

### Measuring example

Mass <i>m</i> in kg	*Force F in N	Spring elongation s in cm
0	0	0
0.3	3	9
0.5	5	15
0.7	7	21
1.0	10	30

<sup>\*</sup>  $F = m \cdot g$  (round values)

### **Evaluation**



At a helical spring, the acting force and the spring elongation are proportional:  $F \sim s$ .

The ratio of the force F and the elongation s is a constant, which is called spring constant D:  $\frac{F}{s} = D$ .