

## Electromagnetism and induction

## Transformers

## Model of a high-current transformer

## Induction furnace

## Objects of the experiment

1. Demonstrating the operation of an induction furnace
2. Measuring the primary current and estimating the secondary current

## Setup



## Evaluation

In the high-current transformer set up for this experiment, the ratio of the numbers of turns of the primary and secondary coil is  $\frac{500}{1}$ .

Due to the current transformation ratio of the high-current transformer ( $\frac{N_1}{N_2} = \frac{I_2}{I_1}$ ), the secondary current is 500 times the primary current.

The large secondary current flowing through the melting ladle leads to a strong warming of the melting ring.

## Apparatus

1 coil with 500 turns, main .....	562 21
1 ring-shaped melting ladle.....	562 20
1 melting ring .....	562 32
1 U-core with yoke .....	562 11
1 clamping device .....	562 12
1 demo-multimeter, passive.....	531 905
2 safety connecting leads, 100 cm, black .....	500 644
1 measuring junction box.....	502 05

## Carrying out the experiment

- Put the melting ring into the melting ladle.
- Switch the coil with 500 turns on, read the primary current from the demo-multimeter, and observe the melting ring.
- When the ring has melted, switch the coil with 500 turns off.

## Observation

After the coil with 500 turns has been switched on, the melting ring starts to melt.

## Measuring example

$N_1$	$N_2$	$I_1$ in A	* $I_2$ in A
500	1	1.3	650

$$* I_2 = \frac{N_1 \cdot I_1}{N_2}$$