

## Determining the electric work of an immersion heater using an AC power meter

### Objects of the Experiment

- Determining the electric work  $W_{el}$  of an immersion heater using an AC power meter
- Determining the electric work  $W_{el}$  measuring the voltage  $U$ , the current  $I$  and the time  $t$
- Determining the change of thermal energy  $E_{th}$

### Principles

Electrical work is the work done on a charged particle by an electric field. In this experiment, electrical work is used to change the thermal energy of water. The supplied electrical work  $W_{el}$  is converted into thermal energy  $E_{th}$  in the heating coil of an immersion heater. This leads to a temperature rise of the water.

The supplied electrical work is determined by:

$$W_{el} = U \cdot I \cdot t \quad (1)$$

$U$ : voltage

$I$ : current

$t$ : time

According to equation (1) the electrical work  $W_{el}$  can be determined by measuring the voltage  $U$ , the current  $I$  and the time  $t$ . In addition, the electrical work can be measured directly by an AC meter (alternating current meter).

The change in thermal energy  $E_{th}$  is determined by:

$$E_{th} = c_W \cdot m \cdot (\vartheta_2 - \vartheta_1) \quad (2)$$

$\vartheta_1$ : temperature at start

$\vartheta_2$ : temperature at end

$c_W$ : specific heat capacity of water

$m$ : mass of water

According to equation (2), measuring the temperature  $\vartheta_1$  at the start and  $\vartheta_2$  at the end of the experiment allows to determine the thermal energy  $E_{th}$ .

The supplied electrical work and the thermal energy can be determined quantitatively in units of wattsecond (Ws) and Joule (J) so that their numerical equivalence can be demonstrated experimentally:  $W_{el} = E_{th}$ .

In the experiment the amount of electrical work  $W_{el}$  is determined that is required to bring 1 liter of water from room temperature to a boil. Three different measuring methods are applied:

- a) by noting the meter reading  $M_1$  of a electricity meter before beginning experiment and the meter reading  $M_2$  when the water boils,
- b) by counting the number of revolutions of the electricity meter disc during the heating process.



Fig. 1: Experimental setup for method c

c) by measuring the voltage  $U$ , the current  $I$  and the time  $t$  during the heating process.

The change in thermal energy  $E_{th}$  is determined from the temperature change between beginning and end of the heating. The resulting values are compared.

### Apparatus

|  |         |
|--|---------|
| 1 Alternating current meter .....          | 560 331 |
| 1 Pair of stand feet .....                 | 301 339 |
| 1 Immersion heater.....                    | 303 25  |
| 1 Plastic beaker, 1000 ml .....            | 590 06  |
| 1 Thermometer, -10 to +110 °C.....         | 382 34  |
| 1 Multimeter LDanalog 20 .....             | 531 120 |
| 1 Multimeter LDanalog 30 .....             | 531 130 |
| 1 Stop-clock I.....                        | 313 07  |
| 4 Safety connecting lead 50 cm black ..... | 500 624 |

**Safety notes**

Use caution when experimenting with dangerous contact voltages:

- Connect the device in circuits only when the mains plug is unplugged from the mains.
- Only connect the device to the AC mains using the power cable
- Connect voltage and current meters to the device using only safety connecting leads (500 420 ff) and only at the safety connection sockets.
- When no ammeter is connected, use only the safety bridging plug supplied with the AC meter to jumper the socket pair to the right.
- Observe the notes contained in the instruction sheet of the alternating current meter (560 331) and the immersion heater (303 25).

**Setup**

- Set up the alternating current meter on the pair of stand feet.
- Jumper the safety socket pair to the right (above the mains socket for the consumer) with the safety bridging plug.
- Connect the alternating current meter to the AC mains (230 V / 50 Hz).
- Pour 1.0 liter water in the plastic beaker.
- Put the immersion heater in the plastic beaker.

**Note:**

The immersion heater is equipped with an overheating protection. **Never** use immersion heater with less water than the minimum immersion level indicated in the operation instructions! Otherwise the heater will permanently fail to work.

**Carrying out the experiment****Method a)**

- Measure the temperature  $\vartheta_1$  of the water.
- Note meter reading  $M_1$  from the counter value display.
- Connect the immersion heater to mains socket for the consumer of the alternating current meter. The immersion heater starts heating the water.
- Heat water till it boils and withdraw mains plug.
- Measure the temperature  $\vartheta_2$  of the water.
- Note meter reading  $M_2$ .

**Method b)**

- Measure the temperature  $\vartheta_1$  of the water.
- Connect the immersion heater to mains socket for the consumer of the alternating current meter. The immersion heater starts heating the water.
- Count the number  $n$  of revolutions of the meter disc during the heating process.
- Heat water till it boils and withdraw mains plug.
- Measure the temperature  $\vartheta_2$  of the water.

**Method c)**

- Withdraw the main plug of the alternating current meter.
- Connect the multimeter LDanalog20 as voltmeter to the safety sockets on the left side of the alternating current meter.
- Remove the safety bridging plug on the right side and connect the multimeter LDanalog30 as ammeter.

- Connect the alternating current meter to the AC mains (230 V / 50 Hz).
- Measure the temperature  $\vartheta_1$  of the water.
- Connect the immersion heater to the mains socket for the consumer of the alternating current meter and simultaneously start stop watch.
- Read off applied voltage  $U$  and current  $I$ .
- Heat water till it boils.
- Withdraw mains plug and simultaneously stop stop watch.
- Measure the temperature  $\vartheta_2$  of the water.

**Measuring example and evaluation****Method a)**

- Meter reading  $M_1$ : 12.70 kWh
- Meter reading  $M_2$ : 12.78 kWh
- Temperature  $\vartheta_1$ : 32 °C
- Temperature  $\vartheta_2$ : 100 °C
- Specific heat capacity of water:  $c_W = 4.18 \frac{\text{kJ}}{\text{kg} \cdot \text{K}}$
- Mass of 1.0 liter water: 1.0 kg

The electrical work  $W_{\text{el}}$  is calculated from the difference  $M_2 - M_1 = W_{\text{el}} = 0.08 \text{ kWh} = 288 \text{ kWs}$

The change in thermal energy  $E_{\text{th}}$  is determined according equation 2:  $E_{\text{th}} = c_W \cdot m \cdot (\vartheta_2 - \vartheta_1) = 284 \text{ kJ}$

**Method b)**

- Number of revolutions of the meter disc  $n$ : 47
- Meter constant  $k$ :  $600 \frac{\text{revs.}}{\text{kWh}}$
- Temperature  $\vartheta_1$ : 32 °C
- Temperature  $\vartheta_2$ : 100 °C

The electrical work  $W_{\text{el}}$  is calculated from the number  $n$  of revolutions and the meter constant  $k$ :

$$W_{\text{el}} = \frac{n}{k} = 0.078 \text{ kWh} = 282 \text{ kWs}$$

The change in thermal energy  $E_{\text{th}}$  is determined according equation 2:  $E_{\text{th}} = c_W \cdot m \cdot (\vartheta_2 - \vartheta_1) = 284 \text{ kJ}$

**Method c)**

- Voltage  $U$ : 220 V
- Current  $I$ : 4.1 A
- Time  $t$ : 5 min 24 s = 0,09 h
- Temperature  $\vartheta_1$ : 32 °C
- Temperature  $\vartheta_2$ : 100 °C

The electrical work  $W_{\text{el}}$  is calculated according equation 1:

$$W_{\text{el}} = U \cdot I \cdot t = 0.081 \text{ kWh} = 292 \text{ kWs}$$

The change in thermal energy  $E_{\text{th}}$  is determined according equation 2:  $E_{\text{th}} = c_W \cdot m \cdot (\vartheta_2 - \vartheta_1) = 284 \text{ kJ}$

**Comparison of the experimental results:**

Electrical work  $W_{\text{el}}$ :

Method 1:  $W_{\text{el}} = 288 \text{ kWs}$

Method 2:  $W_{\text{el}} = 282 \text{ kWs}$

Method 3:  $W_{\text{el}} = 292 \text{ kWs}$

Change in thermal energy  $E_{\text{th}}$ :

$$E_{\text{th}} = 284 \text{ kJ}$$