

Investigation of pressure-dependent line broadening by means of a Michelson interferometer

Objects of the experiments

- Determination of the coherence length of the green spectral line of an Hg spectral lamp and a high pressure mercury lamp
- Determination of the coherence time and line width of the green spectral line of an Hg spectral lamp and a high pressure mercury lamp

Principles

Coherence is the ability of various waves to create stationary interference effects. A temporally stationary interference structure can only be observed when the phase differences between any partial waves around a fixed point change during the observation time by less than 2π . Then the partial waves are called temporally coherent. The maximum time span Δt , during which the phase differences between all parts of the waves change by a maximum of 2π is called the coherence time.

The coherence time Δt_C is directly associated with the spectral width $\Delta\nu$ of the light source. The following applies:

$\Delta t_C = \frac{1}{\Delta\nu}$. By determination of the coherence time Δt_C the spectral width $\Delta\nu$ or $\Delta\lambda$ of a light source, e.g. of individual spectral lines, is determined from the relationships

$$\Delta\nu = \frac{1}{\Delta t_C} \quad \text{or} \quad \Delta\lambda = \frac{1}{c} \frac{\lambda_0^2}{\Delta t_C} \quad (\text{I})$$

Often the coherence length is used instead of the coherence time. This describes the distance

$$\Delta s_C = \frac{c}{n} \Delta t_C \quad (\text{II})$$

the light travels in a medium with a refractive index n during the coherence time.

Typical coherence lengths are several metres for lasers, several centimetres for spectral lines and a few micrometers for white light sources. They can e.g. be measured by means of the Michelson interferometer (see experiment P5.3.4.4).

Natural line widths for spectral lines in the visible spectral range are approx. $\Delta\lambda \approx 10^{-14}$ m. This corresponds to a coherence time of approx. 10^{-8} s or a coherence length of 30 m. In commercial spectral lamps the temperature and pressure situation will, however, lead to a noticeable broadening of the spectral lines.

The dominant effect of the line broadening in spectral lamps is pressure broadening. If during the light emission, impact with a further atom occurs, this leads to a change in the photon energy and/or the phase of the emitted wave and therefore to a change in the line width.




Fig. 1: Experimental setup

This effect is linearly associated with the gas pressure and leads to additional shifting of the spectral lines.

A further part of line broadening is based on the Doppler effect because the atoms move statistically in space during the emission. This broadening mechanism increases linearly with the translation speed of the atoms and therefore with increasing temperatures \sqrt{T} . In the visible spectral range the Doppler broadening exceeds the natural line width by approximately two orders of magnitude.

In the experiment, the coherence lengths of an Hg spectral lamp and an Hg high pressure lamp are determined and compared. For determination of the coherence lengths, the positions of an adjustable planar mirror are determined for which interference can still just be observed. From the difference in the path lengths the coherence length Δs_C can be directly derived, and from this the coherence time Δt_C and the line width $\Delta\nu$ of the spectral line can be determined. The various results for the spectral lamp and the high pressure mercury lamp are discussed.




Fig. 2: The setup of the Michelson interferometer on the optical bench seen from above

- a spectral lamp Hg 100 / high pressure mercury lamp
- b, d iris diaphragm
- c monochromatic filter, yellow-green
- e, g planar mirror with fine adjustment
- f beam splitter
- h fine adjustment drive
- i reduction gears of the fine adjustment drive
- k translucent screen

Safety notes:

Spectral lamp Hg 100

- Connect spectral lamp (451 062) in housing (451 16) only via the universal choke (451 30) to the mains.
- Between the light opening and the optical element (e.g. diaphragm, lens) a minimum distance of 3 cm has to be adhered to in order to prevent overheating.

High pressure mercury lamp

- Caution! The lamp will heat up to temperatures above 100 °C!
- Only connect to the alternating voltage mains (230 V~) via a choke.
- Do not look into the direct or reflected light beam (UV radiation)
- Protect the lamp against impact, falling etc. (risk of breaking)

Observe the operating instructions!

Setup

Note: Carry out the measurements in a room as completely dark as possible.

Optical components with damaged or dirty surfaces can lead to errors in the interference pattern.

Treat the planar mirror and the beam splitter with great care, store them protected from dust and never touch with bare hands.

A experimental setup applying the cross connector (460 342) is shown at the end of the leaflet.

Apparatus

1 spectral lamp Hg 100	451 062
1 housing for spectral lamps	451 16
1 high pressure mercury lamp.....	451 15
1 E27 socket, multi-way plug.....	451 19
1 universal choke 230 V, 50 Hz.....	451 30
1 optical bench with standard profile, 1 m	460 32
1 optics rider 60/50.....	460 373
7 optics riders 90/50	460 374
2 planar mirrors with fine adjustment.....	473 461
1 fine adjustment drive	473 48
1 cantilever arm.....	460 380
1 beam splitter.....	471 88
2 iris diaphragms	460 26
1 monochromatic filter, yellow-green.....	468 07
1 holder with spring clips	460 22
1 translucent screen	441 53
1 saddle base.....	300 11

The experimental setup is shown in figure 2.

For setting up, the steps described below are required:

Installation on the optical bench:

- Mount the Hg 100 spectral lamp (a) in the optics rider 60/50 at one end of the optical bench.
- Mount the reduction gears for the fine adjustment drive (i) with the magnetic base on the gear base and mount it on the other end of the optical bench.
- Clamp a planar mirror (g) at the top end of the fine adjustment drive (h) and mount it in front of the reduction gears (i) on the bench.
- Carefully clamp the universal coupling in the joint head of the micrometer screw of the fine adjustment drive (h)
- Shift the optics rider with the fine adjustment drive (h) and the height of the gear base of the reduction gears (i) in such a way that the coupling rods are neither fully stretched out nor compressed. Otherwise the measurement might become falsified because of shifting of the fine adjustment drive.

- Keep the angle between the individual elements of the couplings as small as possible (and under no circumstance larger than 45°)
- Mount one of the iris diaphragms (**b**) approx. 5 cm behind the spectral lamp Hg 100, the second iris diaphragm (**d**) approx. 25 cm in front of the planar mirror, in such a way that the centres of the diaphragms are at the same height.
- Connect the universal choke to the spectral lamp and switch on; wait for several minutes for the lamp to warm up.

Adjustment of the planar mirror (g)

After the adjustment of the planar mirror (**g**) it should reflect the light from the spectral lamp back onto its own path. Only in this case shifting of the planar mirror (**g**) by means of the fine adjustment drive will not lead to a beam movement.

- Close the iris diaphragms (**b**, **d**) as much as possible.
- Align the planar mirror (**g**) by adjusting the screws on the back so that the reflected beam hits both iris diaphragm (**d**) and iris diaphragm (**b**).

After this the adjustment screws on the planar mirror (**g**) must not be touched again! A change to the adjustment of the planar mirror (**g**) would mean that any shifting of the planar mirror (**g**) with the fine adjustment drive would result in a movement of the reflected beam and therefore the partial beams would no longer superimpose properly.

Beam splitter (f) and the second planar mirror (e)

- Place the beam splitter (**f**) as close as possible to the front of the planar mirror (**g**) so that the complete travel of the fine adjustment drive can be utilised. The mirror side of the beam splitter points in the direction of the spectral lamp.
- Turn the beam splitter (**f**) in such a way that the beam reflected from the planar mirror (**g**) is deflected by 90°.
- Fix the translucent screen (**k**) in the base and place it next to the optical bench so that it is hit at its centre.
- Then clamp the second planar mirror (**e**) in the cantilever arm and mount it on the optical bench in such a way that the planar mirror is hit at its centre by the partial beam reflected from the beam splitter. The distances the planar mirror (**e**) and planar mirror (**g**) are from the beam splitter (**f**) should be approximately equal. On account of the relatively small coherence length of the spectral lamp of a few millimetres, the optical wavelength of the two interferometer arms may be only very slightly different. If necessary adjust the positions of the planar mirror (**e**) and/or the beam splitter (**f**).

Superimposition of the partial beams

For the superimposition adjust only the planar mirror (**e**) (on the cantilever arm)!

- Align the planar mirror (**e**) by adjusting the screws on the back so that the beam is nearly reflected along its own path and after transmission through the beam splitter combines with the first partial beam.
- By fine adjustment of the planar mirror (**e**) the beams of the two interferometer arms can be fully superimposed. For doing this it is useful to cover the beam reflected from the planar mirror (**e**) directly in front of the planar mirror (**e**) partially with a piece of firm paper (e.g. a calling card). The position of the partial beam allowed to

pass can now easily be compared to the position of the beam reflected by the planar mirror (**g**).

- Fully open iris diaphragm (**d**).

Now interference lines should become visible on the screen.

- Carefully adjust the planar mirror (**e**) in such a way that on the screen a system of concentric rings becomes visible in the centre of the lit area.
- Open the iris diaphragm (**b**) sufficiently that the contrast in the interference pattern is not affected.

Carrying out the experiment

During the experiment:

- Avoid mechanical vibration of the optical bench (e.g. do not wobble the table)
- Avoid air flow through the setup (because of flow marks) e.g. from draughts.
- Mark a location on the transparent screen (**k**) where the drifting interference lines can be counted.
- Because of the play in the gearbox, adjust the gearbox knob slowly and uniformly by gently placing the finger onto the lever of the reduction gears (**i**) and in this way until, if necessary with more turns, the interference lines start moving.
- Then give the gearbox knob at least one further turn before starting with the measurement.

Note: if the planar mirror and therefore the interference pattern moves jerkily, the slide bush of the fine adjustment drive needs to be lubricated.

a) Determination of the coherence length of the green Hg spectral line

- Clamp the monochrome filter, yellow-green, in the holder with spring clips and mount it behind the iris diaphragm (**b**) on the optical bench
- Rotate the gearbox knob until the interference pattern is only just visible. Read the setting from the micrometer screw of the fine adjustment drive.
- Rotate the gearbox knob in the other direction to make the interference pattern reappear. Continue rotating until the interference pattern is again only just visible and read the setting from the micrometer screw of the fine adjustment screw.
- Place the planar mirror (**g**) into the so-called white light position e.g. into the centre between the two positions where interference is just visible.

b) Determination of the coherence length of the green Hg spectral line for the high pressure mercury lamp

Caution!

For dismantling allow the spectral lamp to cool down or only touch it with a cloth.

Do not change the setting and position of the planar mirror and the beam splitter!

- Replace the Hg spectral lamp by the high pressure mercury lamp.
- Switch on the universal choke and wait for a few minutes to allow the lamp to heat up.
- Now the interference pattern should be visible again. If necessary check the position of the planar mirror (**e**).

- Open the iris diaphragm (b) sufficiently that the contrast in the interference pattern is just not affected.
- Rotate the gearbox knob until the interference pattern is only just visible. Read the setting from the micrometer screw of the fine adjustment drive.
- Rotate the gearbox knob in the other direction to make the interference pattern reappear. Continue rotating until the interference pattern is again only just visible and read the setting from the micrometer screw of the fine adjustment screw.

Measuring examples and evaluation

a) Determination of the coherence length of the green Hg spectral line for the Hg spectral lamp

The coherence length Δs_C corresponds to the difference between the two optical path lengths and therefore the two mirror positions x_1 and x_2 where interference can still just be observed.

$$\Delta s_C = 2 (x_2 - x_1)$$

The experiment gives $x_1 = 12.5$ mm and $x_2 = 18.7$ mm and therefore the coherence length of $\Delta s_C = 12.4$ mm




Fig. 3: Interference pattern on the screen without filter

Calculation of the coherence time and line width

The coherence time results from formula (II)

$$\Delta t_C = \frac{n}{c} \Delta s_C.$$

With $\Delta s_C = 12.4$ mm one obtains with $n = 1.0$ and $c = 3.0 \cdot 10^8 \frac{\text{m}}{\text{s}}$ for the coherence time $\Delta t_C = 4.1 \cdot 10^{-11}$ s.

The line width of the green Hg spectral line is found from the formula (I): $\Delta \nu = \frac{1}{\Delta t_C}$ or $\Delta \lambda = \frac{1}{c} \frac{\lambda_0^2}{\Delta t_C}$.

For $\Delta t_C = 4.1 \cdot 10^{-11}$ s and $\lambda = 546$ nm this gives $\Delta \nu = 2.4 \cdot 10^{10} \frac{1}{\text{s}}$ or $\Delta \lambda = 2.4 \cdot 10^{-11}$ m.

b) Determination of the coherence length of the green Hg spectral line for the high pressure mercury lamp

The experiment gives $x_1 = 14.7$ mm and $x_2 = 16.5$ mm and therefore the coherence length

$$\Delta s_C = 2 (x_2 - x_1) = 3.6 \text{ mm.}$$

Calculation of the coherence time and line width

With $\Delta s_C = 3.6$ mm one obtains with $n = 1.0$ and $c = 3.0 \cdot 10^8 \frac{\text{m}}{\text{s}}$ for the coherence time $\Delta t_C = 1.2 \cdot 10^{-11}$ s.

The line width of the green Hg spectral line is determined from $\Delta t_C = 1.2 \cdot 10^{-11}$ s and $\lambda = 546$ nm:

$$\Delta \nu = 8.3 \cdot 10^{10} \frac{1}{\text{s}} \text{ or } \Delta \lambda = 8.3 \cdot 10^{-11} \text{ m.}$$

As described under Principles, the coherence lengths for the Hg spectral lamp and the high pressure Hg lamp have values of only a few millimetres. In addition they differ by a factor of 3.4 which is mainly caused by the different pressures in the lamps.

Note:

If the monochromatic filter, yellow-green, is removed from the path of the beam and the planar mirror (g) is shifted carefully by rotating the gearbox knob, it is apparent that in certain regions interference images are obtained with a strong light-dark contrast while in other positions a coloured ring system is obtained.

In the case of double lines, the superposition in the Michelson interferometer leads to further variations in intensity. This effect can be utilised in order to determine the line separation in the case of double lines (see experiment P5.3.4.6).

Alternative Setup with cross connector (460 342)




Fig. 4: The setup of the Michelson interferometer on the optical bench seen from above

- a** spectral lamp Hg 100 /
high pressure mercury lamp
- b, d** iris diaphragm
- c** monochromatic filter, yellow-green
- e** planar mirror with fine adjustment in extension rod
- f** beam splitter
- g** planar mirror with fine adjustment
- h** fine adjustment drive
- i** reduction gears of the fine adjustment
- k** translucent screen

Connect the optical benches with the cross connector as shown in Fig. 4. Then, the setup of the optical components is similar to the setup on a single optical bench.

Apparatus

1 spectral lamp Hg 100	451 062
1 housing for spectral lamps	451 16
1 high pressure mercury lamp.....	451 15
1 E27 socket, multi-way plug.....	451 19
1 universal choke 230 V, 50 Hz.....	451 30
3 optical bench with standard profile, 0.5 m	460 335
1 cross connector.....	460 342
1 optics rider 60/50.....	460 373
7 optics riders 90/50.....	460 374
1 extension rod.....	460 385
2 planar mirrors with fine adjustment	473 461
1 fine adjustment drive	473 48
1 beam splitter.....	471 88
2 iris diaphragms.....	460 26
1 monochromatic filter, yellow-green.....	468 07
1 holder with spring clips.....	460 22
1 translucent screen.....	441 53
1 saddle base.....	300 11

