

## Malus' Law

### Objectives of the Experiment

- Understanding the polarization of light.
- Investigations of the functional principle of a polarizer and an analyzer.
- Qualitative validation of Malus' Law.

### Principles

The phenomena of electrodynamics and magnetism were unified by the Maxwell equations for *electromagnetism*. The Maxwell equations in a vacuum can be solved by wave functions and thus describe electromagnetic waves. Such waves are transverse waves, i.e. the propagation direction  $\vec{k}$  of the wave is always perpendicular to the electric field  $\vec{E}$  and the magnetic flux density  $\vec{B}$ :

$$\vec{k} \cdot \vec{E} = \vec{k} \cdot \vec{B} = 0. \quad (1)$$

This property of the electromagnetic waves is instrumental for the phenomenon of polarization. Light is an electromagnetic wave as well. Polarization of light is the direction of the vector  $\vec{E}$  or  $\vec{B}$  in the plane perpendicular to  $\vec{k}$ . The unpolarized light is composed of superimposed waves with all possible directions of  $\vec{E}$  and  $\vec{B}$  combined. If the direction is constant in time, we have *linear polarization*. The  $E$  field periodically changes its sign and absolute value. For *circular polarization*, the vector  $\vec{E}$  rotates in the plane perpendicular to the direction of propagation  $\vec{k}$  and has a constant absolute value.

The thermally generated light (lightbulb) is unpolarized and consists of light waves that are linearly polarized in all directions in the plane perpendicular to  $\vec{k}$ . When passing through a polarization filter with polarization direction  $\vec{v}$ , each vector  $\vec{E}$  is reduced to its projection

$$E_v = \vec{E} \cdot \vec{v} = E \cos \theta \quad (2)$$

in the light beam. Here,  $\theta$  is the angle between  $\vec{E}$  and  $\vec{v}$ . The light intensity averaged over all angles  $\theta$  is halved. The light is also completely linearly polarized, i.e. the entire  $E$  field

$$\vec{E}_v = E_v \vec{v} \quad (3)$$

is parallel to  $\vec{v}$ . The intensity

$$I_0 = (E_v)^2 \quad (4)$$

Of the polarized light is a square functions of the  $E$  field.




Figure 1: A light beam is linearly polarized in the polarization direction upon passing through a polarization filter.

When passing through a further polarization filter with polarization direction  $\vec{w}$ , forming an angle  $\alpha$  with  $\vec{v}$ , the vector  $\vec{E}_v$  is reduced to its projection

$$E_w = \vec{E}_v \cdot \vec{w} = E_v \cos \alpha \quad (5)$$

The intensity

$$I = (E_w)^2 = (E_v)^2 \cos^2 \alpha \quad (6)$$

is again a square function of the  $E$  field. Figure 1 schematically shows the functional principle of a polarization filter. Substituting equation (4) into equation (6) yields Malus' Law follows

$$I = I_0 \cos^2 \alpha. \quad (7)$$

In the present experiment, Malus' Law will be confirmed quantitatively.

**Apparatus**

1 Halogen lamp, 12 V, 50/100 W.....	450 64
1 Halogen bulb, 12 V/100 W, G6.35.....	450 63
1 Picture slider.....	450 66
1 Transformer, 2 ... 12 V, 120 W.....	521 25
1 Iris diaphragm.....	460 26
2 Polarization filters.....	472 401
1 Lens in frame, f=100 mm.....	460 03
1 Solar cell, STE 2/19.....	578 62
1 Holder for plug-in elements.....	460 21
1 Digital multimeter 3340.....	531 183
1 Small optical bench.....	460 43
6 Leybold multiclamp.....	301 01
2 Stand base, V-shaped, large.....	300 01
2 Connecting lead, 32 A, 100 cm, black.....	501 33
1 Connecting lead 19A,	
3 100 cm, red/blue, pair.....	501 46

**Safety Note**

Watch for heat development of the halogen lamp and do not cover the ventilation slots!

**Experiment set-up**

The experiment set-up is shown in Figure 2.

- Attach the small optical bench to the stand base.
- Carefully align the optical bench in a horizontal position.
- Insert the halogen bulb into halogen lamp and mount the latter to one end of the optical bench. Use Leybold multiclamps for all optical devices.
- Mount the picture slider with the heat protection filter for absorption of the infrared component in the halogen light to the halogen lamp.
- Mount the iris diaphragm on the optical bench in front of the halogen lamp so that the distance between the picture slider and the iris diaphragm is approximately 2 - 5 mm.
- Adjust the aperture of the iris diaphragm to a diameter of 2 - 4 cm.
- Connect the halogen lamp to the transformer with the black connecting leads.
- Mount the holder for the plug elements at the other end of the optical bench and insert the solar cell.
- Connect the solar cell to the digital multimeter to measure the current and set it to the finest current measurement.
- Move the Leybold multiclamp with the lens  $f = +100$  along the optical bench until the image of the iris diaphragm is focussed on the solar cell. If necessary, hold up a piece of white paper.
- Mount the two polarizing filters between the lens and the solar cell and set them to  $0^\circ$ .



Figure 3: Experiment set-up, Malus' Law.

**Performing the experiment**

In the experiment, the current of the solar cell is measured, which is proportional to the light intensity. Equation (7) thus also applies to the measured currents.

- Set both polarization filters to  $0^\circ$  and record the solar cell current. Thus,  $\alpha = 0^\circ$  and the current value corresponds to  $I_0$ .
- Adjust the polarization filter on the solar cell from  $0^\circ$  to  $180^\circ$  in increments of  $10^\circ$ . Record the angle difference  $\alpha$  as well as the solar cell current  $I$ .
- The current value at  $\alpha = 90^\circ$  corresponds to the background and must be subtracted from each current value for the analysis.

**Measurement example**

For the measurement, the angle difference  $\alpha$  between the polarization filters was varied from  $0^\circ$  to  $180^\circ$  in increments of  $10^\circ$ . The measurement example is shown in Figure 3.




Figure 2: Measurement example, Malus' Law: (□) measured values, (—) equation (7) for measured current  $I_0 = 75 \mu\text{A}$ .