

Rotation of the plane of polarization with quartz

Objects of the experiments

- Observing the rotation of the plane of polarization by quartz in an arrangement of two crossed polarizers.
- Distinguishing between right-handed, left-handed and non-rotating (i. e. cut parallel to the crystal optic axis) quartz.




Fig. 1 Experimental setup for observing the rotation of the plane of polarization by quartz.

- a Halogen lamp housing
- b Light filter (in picture slider)
- c Polarizer
- d Quartz crystal (in holder with spring clips)
- e Analyzer
- f Lens
- g Observing screen

Principles

If a plane-parallel quartz plate which is cut perpendicularly to a crystal direction called the optic axis is put into the ray path between two crossed polarizers, the field of view is lit up (*F. Arago*, 1811). For monochromatic light darkness is re-established if the analyzer is rotated by a certain angle α . That means the quartz plate has rotated the plane of polarization of the light by the angle α .

The rotatory power of the quartz crystal is due to its helical crystal structure, which leads to different phase velocities for the propagation of right-circularly and left-circularly polarized light in the crystal. Linearly polarized light which enters the crystal can be decomposed into a right-circularly and a left-circularly polarized partial wave. The two partial waves propagate at different phase velocities so that a phase difference arises, which is proportional to the distance covered in the crystal. After the two partial waves have covered this distance, their superposition results in a linearly polarized wave whose direction of polarization is rotated relative to the original wave. The angle of rotation α thus is proportional to the distance covered in the crystal, i. e. to the thickness d of the crystal.

There are right-handed and left-handed quartzes. In a right rotation the plane of polarization is rotated clockwise as seen by an observer looking in the direction opposite to the light ray propagation, in a left rotation the plane of polarization is rotated anticlockwise. As the angle of rotation strongly depends on the wavelength of the light, single-coloured light is used in the experiment.

Apparatus

1 quartz, parallel	472 62
1 quartz, right-handed	472 64
1 quartz, left-handed	472 65
1 halogen lamp, 12 V/100 W	450 63
1 halogen lamp housing 12 V, 50/100 W	450 64
1 picture slider for halogen lamp housing	450 66
1 mercury light filter, yellow	468 30
1 transformer 2 ... 12 V	521 25
2 polarization filters	472 401
1 lens, $f = + 100$ mm	460 03
1 holder with spring clips	460 22
1 translucent screen	441 53
1 small optical bench	460 43
1 stand base, V-shape, 28 cm	30001
6 Leybold multiclips	301 01
Connecting leads with 2.5 mm ² cross section	

Evaluation and results

Viewed opposite to the direction of the light ray, the right-handed quartz crystal rotates the plane of polarization by 32.5° to the right whereas the left-handed quartz crystal rotates it to the left.

Setup

The experimental setup is illustrated in Fig. 1.

- Mount the components on the small optical bench according to Fig. 1, where the position of the left edge of the Leybold multiclips is given.
- Align the polarization filters so that their scales point towards the observing screen, and set them both to 90° .
- Set up the halogen lamp housing for 100 W operation (use the reflector, see instruction sheet for halogen lamp housing), and insert the mercury light filter in the picture slider in front of the exit aperture.
- Align the halogen lamp with the adjusting rod (**a1**) of the lamp housing and shift the lens on the optical bench so that the field of view on the observing screen is uniformly illuminated.
- Set the analyzer to 0° (maximum darkness of the field of view).

Carrying out the experiment

- Insert the parallel quartz in the holder with spring clips so that it is centred, and adjust maximum darkness with the analyzer.
- Insert the right-handed quartz in the holder with spring clips so that it is centred, and adjust maximum darkness with the analyzer.
- Insert the left-handed quartz in the holder with spring clips so that it is centred, and adjust maximum darkness with the analyzer.

Measuring example

Table 1: Angular position of the analyzer for maximum darkness of the field of view (polarizer: 90°)

object	thickness	angular position
without quartz		0°
quartz, parallel		0°
quartz, right-handed	1,5 mm	$+32,5^\circ$
	4 mm	$+87^\circ$
quartz, left-handed	1,5 mm	$-32,5^\circ$
	4 mm	-87°